

PLR compressed air lab stirrers

Accessories for compressed air lab stirrers

Bossheads

Bosshead KR 260

Clamping width max. 36 mm


Stainless steel bosshead KR 360

Bosshead with prism clamping on both sides, clamping range 12 to 36 mm, heavy design all in stainless steel.

Bosshead KK-16, swivelling

Clamping width 16 mm

Type	Code-No.
KR 260	40 26446 00386 9
KR 360	40 26446 00815 4
KK-16	40 26446 00325 8



KR 260



KR 360



KK-16



Clamp holder SP

With bosshead KR 260, strap length 150 mm

Type	Description	Code-No.
SP	Clamp holder with bosshead KR 260	40 26446 00792 8



Stirring shaft protection RWS

Prevents injuries from rotating stirrers during use. Length adjustable 190-310 mm, incl. support arm L=275 mm. Bosshead not included.

Type	Description	Code-No.
RWS	Stirring shaft protection	40 26446 00095 0



Service unit type WE-2018

The service unit with filter, regulator and lubricator function provides both clean air with the recommended air purity class at the output and lubrication of motors.

Service unit combination consisting of filter regulator with pressure gauge, rotary knob with detent, mist lubricator. Pneumatic connection G1/4".

Type	Description	Code-No.
WE-2018	Service unit	40 26446 00751 5