

PLR compressed air lab stirrers

Accessories for compressed air lab stirrers

Finger-operated oiler OT

For manual lubrication of compressed air coming into the motor.

Type	Description	Code-No.
OT	Finger-operated oiler	40 26446 00047 9



Motor oil

Cylinder oil

Resin-free and non-corrosive oil for lubrication of compressed air stirrers.

Motor oil for the food industry LT-Oil

Fully synthetic food grade oil according to NSF H1, FDA 21 CFR 178.3570, LMBG §31 para. 1. Operating temperature -60°C to +135°C.

Type	Description	Code-No.
Oil	Cylinder oil 500 ml	40 26446 00081 3
LT-Oil	Food grade oil 500 ml	40 26446 00524 5



Speed regulation valve V

Fine control valve for stepless speed regulation.

Type	Description	Code-No.
V	Speed regulation valve	40 26446 00046 2

