

Pharma/Food grade stirrers

Accessories for pharma/food grade stirrers



Connection couplings VK

Fixed-diameter connection between motor and stirring shafts. Suitable for use in potentially explosive atmospheres (ATEX).

Type	Description	Material no.	Receiver Ø mm	Code-No.
VK 12 x 10	Connection couplings VK	1.4404	12 x 10 mm	40 26446 00525 2
VK 12 x 12	Connection couplings VK	1.4404	12 x 12 mm	40 26446 00526 9
VK 12 x 14	Connection couplings VK	1.4404	12 x 14 mm	40 26446 00647 1

Stands

Floor stand BS

With angled foot, powder-coated or stainless steel, side length 600 mm, stainless steel rod Ø 34 mm

Wall bracket WS

Stainless steel rod Ø 34 mm, clearance from wall: 200 mm

Type	Description	Length	Code-No.
BS 1	Plastic coating / stainless steel	1000 mm	40 26446 00086 8
BS 2	Plastic coating / stainless steel	1500 mm	40 26446 00083 7
BS 3	All in stainless steel	1000 mm	40 26446 00826 0
WS	Stainless steel rod	1000 mm	40 26446 00087 5



BS1 and BS2



BS3



WS

Bossheads


Bosshead KR 260

Clamping width max. 36 mm

Stainless steel bosshead KR 360

Bosshead with prism clamping on both sides, clamping range 12 to 36 mm, heavy design all in stainless steel.

Type	Code-No.
KR 260	40 26446 00386 9
KR 360	40 26446 00815 4



KR 260



KR 360